

## Pansw EKG cable

**DIN 3.0 connector**

**banana 4.0 end**

**Snap end**

**Clips end**



Note: (within this catalogue)

1, integrated means complete cable from machine to electrode; trunk cable means cable from machine to the yoke; leadwire means the wires part from yoke to the electrodes;

2, EKG cable: means 10 leads patient cable for ECG machine or defibrillator;

### Schiller

Part No	OEM No	Description	Picture	Compatibility
EKG01IOD		integrated; IEC, DB15M-Din3.0		Schiller ECG: AT-102, CS-200, Cardiovit AT-10 plus, AT-104 PC, AT-110
EKG02AOD		integrated; AHA, DB15M-Din3.0		Schiller: AT1 ,AT2, AT3, AT5, AT6, AT10, AT60,AT104 PC,CS100,CS6, CS 6/12,AT4, AT104, AT101, AT102; Avionics(Del Mar): 910 / 920 / 930/ 940 / 942 / 950 / 960 / 970 / 980
EKG03IOB	AT-1	integrated; IEC,DB15M-Banana4.0		Custo-med: Custo-norm Cardioline: E1, E1R, E;Delta 3; Delta 300; C3 Cardiorapid: K11, K111, K112, K131, K31, K33, K36 Cardiotek: EP Tracer Dr Lee: ECG-310A; ECG-310B; Del mar: 910,920,930,940,942,950,

EKG04AOB		integrated; AHA,DB15M-Banana4.0	960,970,980 Elettronica Trentina/Carditte: Autoruler 12/1; 12/3; Microruler 12/1; Microruler 12/3; Microtel: Excel 101; Excel 103; Excel 106;
EKG05IOC		integrated; IEC,DB15M-Clip end	Ergoline: Autoruler, Autoscript 6/12 Cardiette, EK 3003 / 3012, Ergoline Esaote Biomedica: EKG P120 Fidelity: LP 3000 Innomed: Heart Mirror1, Heart Mirror3; Heart Mirror3D; Heart Screen 80G; Heart Screen 112C; Heart Screen 112D;
EKG06AOC		integrated; AHA,DB15M-Clip end	Kontron: EKG P80 Progetti: EPG Smarty; EPG Pocket; EPG Junior; EPG View; EPG 6 View; EPG Plus; EPG 612; EPG 13
EKG07IOS		integrated; IEC,DB15M-snap end	Reynolds: Cardiotrack 12 Schiller:AT3, AT5, AT6,AT10, AT60, AT104 PC, AT104, CS6, CS12,CS100, CS200 Von Berg: Bioset 9000
EKG08AOS		integrated; AHA,DB15M-snap end	
EKG09ITT		ecg trunk cable, 10 leads, IEC; use with EKG-10010D. Etc	
EKG09ATT		ecg trunk cable, 10 leads, AHA color; use with EKG-10010L, etc.	







## HP

Part No	OEM No	Description	Connector	Compatibility
EKG10IOD		integrated; IEC,DB15M-Din3.0		
EKG11AOD		integrated; AHA,DB15M-Din3.0		
EKG12IOB	M3703 A	integrated; IEC,DB15M-Banana4.0		
EKG13IOC		integrated; IEC,DB15M-clips		




M1770A(Pagewriter 200i-300i);  
M171A (Pagewriter 200); M1772A  
(Pagewriter 100);

EKG14AOS	M3702 A	integrated;AHA, DB15M-Banana4.0		
EKG15ILC		HP trim 10 lead set, clip end,IEC		
EKG16ALC		HP Pagewriter trim 10 lead set, clip end, AHA/AAMI		<i>Compatible with: PageWriter Touch, 860284, PageWriter Trim Cardiographs, 860286, 860288, 860290, 860297</i>
EKG17ALB	98980 31291 61(AM MI/AH A)	HP Pagewriter trim 10 lead set, banana 4.0 end, AHA/AAMI		
EKG18ILB	98980 31291 91 (IEC)	HP Pagewriter trim 10 lead set, banana 4.0 end; IEC		
EKG19ILB		HP Pagewriter TRIM limb 4 lead set, banana 4.0 end; IEC		IEC 989803129171
EKG20LCHE		HP Pagewriter trim Chest 6 lead set, banana 4.0 end; IEC		IEC 989803129181
EKG20LLIM		HP Pagewriter TRIM limb 4 lead set, banana 4.0 end; AHA/AAMI		AAMI 989803129141
EKG21ALB		HP Pagewriter trim Chest 6 lead set, banana 4.0 end; AHA/AAMI		AAMI 989803129151
EKG22DATA		USB data wire for HP trim ECG acquisition box; original part number: 453564034571		<i>compatible with: Pagewriter Touch (860284) - Pagewriter Trim I , II , III (860286- 860287-860288- 860290)</i>

EKG23ATT		HP DB15 pin connector trunk cable, AHA		
EKG23ITT		HP DB15 pin connector trunk cable, IEC		
EKG24AOC		integrated; AHA; DB15M-clips		Philips Medical: M1770 M1771/M1772
EKG25IOS		integrated; IEC; DB15M-snap		
EKG26AOS		integrated; AHA, DB15M, snap		

### Kenz

Part No	OEM No	Description	Connector	Compatibility
EKG27IOD	PC-104	integrated; IEC, DDK16P-Din3.0		Kenz Delta plus, delta 3plus, delta 60 plus, Cardioline Delta 1- Delta 1 Plus Delta 3-Delta 3 Plus Delta 30 Delta 60-Delta 60 Plus
EKG28AOD		integrated; AHA, DDK16P-Din3.0		
EKG29IOB		integrated; IEC, DDK16P-Banana4.0		
EKG30AOB		integrated; AHA, DDK16P-Banana4.0		
EKG31IOS		integrated; IEC, DDK16P-snap		
EKG32AOS		integrated; AHA, DDK16P-snap		
EKG33IOC		integrated; IEC, DDK16P-clip		
EKG34AOC		integrated; AHA, DDK16P-clip		


### Cardioline; Kenz

Part No	OEM No	Description	Connector	Compatibility
EKG35IOD		integrated; IEC, DB15M-Din3.0		Cardioline 108/110,1203,1205; Suzuken kenz
EKG36AOD		integrated; AHA, DB15M-Din3.0		
EKG37IOB	PC-109	integrated; IEC, DB15M-Banana4.0		
EKG38AOB		integrated; AHA, DB15M-Banana4.0		
EKG39IOS		integrated; IEC, DB15M-snap		
EKG40AOS		integrated; AHA, DB15M-snap		
EKG41IOC		integrated; IEC, DB15M-clip		
EKG42AOC		integrated; AHA, DB15M-clip		



### Nihon Kohden

Part No	OEM No	Description	Connector	Compatibility
EKG43IOD		integrated; IEC, DB15M-Din3.0		Nihon Kohden ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, DongJiang ECG-11A, ECG-11B, ECG-11D, ECG-32A, (II), ECG-120
EKG44AOD		integrated; AHA, DB15M-Din3.0		
EKG45IOB		integrated; IEC, DB15M-Banana4.0		
EKG46AOB		integrated; AHA, DB15M-Banana4.0		
EKG47IOC		integrated; IEC, 10 leads clip end		
EKG48AOC		integrated; AHA, DB15M-clip		
EKG49IOS		integrated; IEC, DB15, snap		
EKG50AOS		integrated; AHA, DB15M-snap		


### Shanghai Kohden

Part No	OEM No	Description	Connector	Compatibility
EKG51IOD		integrated; IEC, DB15M-Din3.0		Shanghai Kohden: ECG-6511, ECG-6151, 8110P/K, 6353
EKG52AOD		integrated; AHA, DB15M-Din3.0		
EKG53IOB		integrated; IEC, DB15M-Banana4.0		
EKG54AOB		integrated; AHA, DB15M-Banana4.0		
EKG55IOC		integrated; IEC, DB15M-snap		
EKG56AOC		integrated; AHA, DB15M-snap		
EKG57IOS		integrated; IEC, DB15M-clip		
EKG58AOS		integrated; AHA, DB15M-clip		

## Burdick


Part No	OEM No	Description	Connector	Compatibility
EKG59IOD	EK-10	integrated; IEC, DB15M-Din3.0		Burdick EK-10
EKG60AOD		integrated; AHA, DB15M-Din3.0		
EKG61IOB		integrated; IEC, DB15M-Banana4.0		
EKG62AOB		integrated; AHA, DB15M-Banana4.0		
EKG63IOS		integrated; IEC, DB15M-snap		
EKG64AOS		integrated; AHA, DB15M-snap		
EKG65IOC		integrated; IEC, DB15M-clip		
EKG66AOC		integrated; AHA, DB15M-clip		
Part No	OEM No	Description	Connector	Compatibility
EKG67IOD		integrated; IEC, DB15M-Din3.0		Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus, Eclipse LE II and ATRIA PLUS
EKG68AOD		integrated; AHA, DB15M-Din3.0		
EKG69IOB		integrated; IEC, DB15M-Banana4.0		
EKG70AOB		integrated; AHA, DB15M-Banana4.0		
EKG71IOS		integrated; IEC, DB15M-snap		
EKG72AOS		integrated; AHA, DB15M-snap		
EKG73IOC		integrated; IEC, DB15M-clip		
EKG74AOC		integrated; AHA, DB15M-clip		

## Fukuda

Part No	OEM No	Description	Connector	Compatibility
EKG75IOD	FX-101	integrated; IEC, DB15M-Din3.0		Fukuda FCP-7411, FX-7402, FX-7202, FX-2111,FX-2155, FX-3010, FX-7102
EKG76AOD		integrated; AHA, DB15M-Din3.0		
EKG77IOB		integrated; IEC, DB15M-Banana4.0		
EKG78AOB		integrated; AHA, DB15M-Banana4.0		
EKG79IOS		integrated; IEC, DB15M-snap		
EKG80AOS		integrated; AHA, DB15M-snap		
EKG81IOC		integrated; IEC, DB15M-clip		
EKG82AOC		integrated; AHA, DB15M-clip		

EKG83IOD		integrated; IEC, DB15M-Din3.0		<i>KP-309/KP-405/KP-500/ C100 C110 C120 C300 C310 C550 C600 C600W C610 C610W C625 a 2000 series; a 1000</i>
EKG84AOD		integrated; AHA, DB15M-Din3.0		
EKG85IOB		integrated; IEC, DB15M-Banana4.0		
EKG86AOB		integrated; AHA, DB15M-Banana4.0		
EKG87IOS		integrated; IEC, DB15M-snap		
EKG88AOS		integrated; AHA, DB15M-snap		
EKG89IOC		integrated; IEC, DB15M-clip		
EKG90AOC		integrated; AHA, DB15M-clip		
EKG91IOD		integrated; IEC, DB15M-Din3.0		<i>KP-309D/KP-405D/KP-500D/ C100 C110 C120 C300 C310 C550 C600 C600W C610 C610W C625 a2000 series; a 1000; (EKG-10009Q---EKG-10009x are with 20k built-in resistor for defibrillation protection)</i>
EKG92AOD		integrated; AHA, DB15M-Din3.0		
EKG93IOB		integrated; IEC, DB15M-Banana4.0		
EKG94AOB		integrated; AHA, DB15M-Banana4.0		
EKG95IOS		integrated; IEC, DB15M-snap		
EKG96AOS		integrated; AHA, DB15M-snap		
EKG97IOC		integrated; IEC, DB15M-clip		
EKG98AOC		integrated; AHA, DB15M-clip		

### GE-Marquette


Part No	OEM No	Description	Connector	Compatibility
EKG99IOB	20298 90-001	Integrated ECG cable, 10 leads, banana end, IEC color coding		<i>Hellige: MicroSmart, MAC 500/1100 ; MAC 1200 / 1200 ST</i>
EKG100AOB	20298 93-001	Integrated ECG cable, 10 leads, banana end, AHA/AAMI color coding		
EKG101IOC		Integrated ECG cable, 10 leads, clips end, IEC color coding		
EKG102AOC		Integrated ECG cable, 10 leads, clips end, AHA/AAMI color coding		
EKG103ITT	22341	Trunk cable, IEC		




EKG104ATT	808	Trunk cable, AHA		
EKG105ILB		Multi-link leadwires, IEC, Banana		
EKG106ALB		Multi-link leadwires, AHA, Banana		
EKG107IOB		integrated; IEC, Banana 4mm		<i>Dash 4000</i>
EKG108IOC		integrated; IEC, clip		
EKG109AOB		integrated; AHA, Banana 4mm		
EKG110AOC		integrated; AHA, clip		
EKG111ITT		Trunk cable, IEC; original part number: (IEC: E9002ZJ, 416035-002)		<i>For GE Eagle, Solar, Dash , Tram MAC-Lab Cath Lab System, Datex-Ohmeda S/5 FM, Original P/N: (AAMI: E9002ZH, 416035-001) (IEC: E9002ZJ, 416035-002)</i>
EKG112ATT		Trunk cable, AHA; original part number: (AAMI: E9002ZH, 416035-001)		
EKG113ILB		10 leads set, IEC, banana, for EKG-10010G		
EKG114ILC		10 leads set, IEC, clip; for EKG-10010G		
EKG115ALB		10 lead set, AHA, banana, for EKG-10010O		
EKG116ALC		10 leads set, AHA, clip; for EKG-10010O		

EKG-10010J	20085 94-004	GE Seer Light one piece, 7 leads, snap		<i>Holter GE Seer Light</i>
------------	-----------------	--	--	-----------------------------

**NEC**

<b>Part No</b>	<b>OEM No</b>	<b>Description</b>	<b>Connector</b>	<b>Compatibility</b>
EKG118IOD		integrated; IEC, NEC15P-Din3.0		



EKG119AOD		integrated; AHA, NEC15P,Din3.0		NEC:K-3321,3330,3370
EKG120IOB		integrated; IEC, NEC15P-Banana4.0		
EKG121AOB		integrated; AHA,15P-Banana4.0		
EKG122IOS		integrated; IEC, NEC15P-snap		
EKG123AOS		integrated; AHA, NEC15P-snap		
EKG124IOC		integrated; IEC, NEC15P-clip		
EKG125AOC		integrated; AHA, NEC15P-clip		

## Suction Electrode

Part No	OEM No	Description	Connector	Patient group
EKG126CHE	21714 401	Suction Electrode use with DIN 3.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;		<i>adult</i>
EKG127CHE		Suction Electrode use with banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;		<i>adult size</i>
EKG128CHE		suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;		<i>adult size</i>
EKG129CHE	20298 92-001; Philips 40418 A	Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 26mm; metal electrode diameter: 21 mm;		<i>adult size</i>
EKG130CHE		Suction Electrode use with DIN 3.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;		<i>pediatric size</i>
EKG131CHE		Suction Electrode use with banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;		<i>pediatric size</i>

EKG132CHE		suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;		<i>pediatric size</i>
EKG133CHE		Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 23mm; metal electrode diameter: 17mm;		<i>pediatric size</i>

### Limb ECG Electrode

Part No	OEM No	Description	Connector	Compatibility
EKG134LIM		Limb Electrode for DIN 3.0 connector; clamps size: 30x130mm;		<i>for adult use</i>
EKG135LIM		Limb Electrode for banana 4.0 connector; clamps size: 30x130mm;		<i>for adult use</i>
EKG136LIM		limb Electrode for DIN 3.0 connector and banana 4.0 connector; clamps size: 30x130mm;		<i>for adult use</i>
EKG137LIM	20298 91-001; Philips 40494 E	Multi-Use Limb Electrode; for snap/DIN/banana/clip connectors; universal use; clamps size: 30x130mm;		<i>for adult use</i>
EKG138LIM		Multi-Use Limb Electrode; for snap/DIN/banana/clip connectors; universal use; clamps size: 20x85mm		<i>for pediatric use</i>

**SIEMENS**

EKG139IOD		integrated; IEC, DB15M-Din3.0		<i>Cardiostat 31</i>
EKG140AOD		integrated; AHA, DB15M-Din3.0		
EKG141IOB		integrated; IEC, DB15M-Banana4.0		
EKG142AOB		integrated; AHA, DB15M-Banana4.0		
EKG143IOS		integrated; IEC, DB15M-snap		
EKG144AOS		integrated; AHA, DB15M-snap		
EKG145IOC		integrated; IEC, DB15M-clip		
EKG146AOC		integrated; AHA, DB15M-clip		

**Bionet**

EKG147IOD		integrated; IEC, DB15M-Din3.0		<i>Cardiocare 2000, Cardiotouch, Medical Econet Cardio M y Cardio M Plus</i>
EKG148AOD		integrated; AHA, DB15M-Din3.0		
EKG149IOB		integrated; IEC, DB15M-Banana4.0		
EKG150AOB		integrated; AHA, DB15M-Banana4.0		
EKG151IOS		integrated; IEC, DB15M-snap		
EKG152AOS		integrated; AHA, DB15M-snap		
EKG153IOC		integrated; IEC, DB15M-clip		
EKG154AOC		integrated; AHA, DB15M-clip		

**MORTARA**

EKG155IOD		integrated; IEC, DB15M-Din3.0		<i>MORTARA INSTRUMENT:ELI 100 MORTARA INSTRUMENT:ELI 100 MORTARA INSTRUMENT:ELI 200 MORTARA INSTRUMENT:ELI 50</i>
EKG156AOD		integrated; AHA, DB15M-Din3.0		
EKG157IOB		integrated; IEC, DB15M-Banana4.0		
EKG158AOB		integrated; AHA, DB15M-Banana4.0		
EKG159IOS		integrated; IEC, DB15M-snap		
EKG160AOS		integrated; AHA, DB15M-snap		
EKG161IOC		integrated; IEC, DB15M-clip		
EKG-10018H	9293-022-50;	integrated; AHA, DB15M-clip		
EKG162IOD		integrated; IEC, DB15M-Din3.0		
EKG163AOD		integrated; AHA, DB15M-Din3.0		

EKG164IOB		integrated; IEC,DB15M-Banana4.0	<i>MORTARA INSTRUMENT: ELI 10/ELI150/ELI250/ ELI210/ELI350</i>
EKG165AOB		integrated; AHA,DB15M-Banana4.0	
EKG166IOS		integrated; IEC, DB15M-snap	
EKG167AOS		integrated; AHA, DB15M-snap	
EKG168IOC		integrated; IEC, DB15M-clip	
EKG169AOC		integrated; AHA, DB15M-clip	

### Cardioline; Cardiette

Part No	OEM No	Description	Connector	Compatibility
EKG170IOD		integrated; IEC, DB15M-Din3.0		<i>AR-600, AR-1200, part number 63050025</i>
EKG171AOD		integrated; AHA,DB15M-Din3.0		
EKG172IOB		integrated; IEC,DB15M-Banana4.0		
EKG173AOB		integrated; AHA,DB15M-Banana4.0		
EKG174IOS		integrated; IEC, DB15M-snap		
EKG175AOS		integrated; AHA, DB15M-snap		
EKG176IOC		integrated; IEC, DB15M-clip		
EKG177AOC		integrated; AHA, DB15M-clip		

### Edan

Part No	OEM No	Description	Connector	Compatibility
EKG178IOD		integrated; IEC, DB15M-Din3.0		<i>SE-12 Express SE-3 SE-601A</i>
EKG179AOD		integrated; AHA,DB15M-Din3.0		
EKG180IOB		integrated; IEC,DB15M-Banana4.0		
EKG181AOB		integrated; AHA,DB15M-Banana4.0		
EKG182IOS		integrated; IEC, DB15M-snap		
EKG183AOS		integrated; AHA, DB15M-snap		
EKG184IOC		integrated; IEC, DB15M-clip		
EKG185IOC		integrated; AHA, DB15M-clip		

### Hellige

Part No	OEM No	Description	Picture	Comaptible
---------	--------	-------------	---------	------------

EKG186IOD		integrated; IEC, DB15M-Din3.0		Bosch:603D,603ER,603S,606 Cambridge: CM 6000, CM 6400 Custo-med: Custocard DEGO: Cardiomed 3 Hellige: EK36,EK41,EK43,EK53, EK 56, EK 403, EK 413, EK 512, EK 512E,EK 512 P, Cardiosmart, CardioSys, Corna Siemens: Cardiostat1/3T/703/710,Cardina C, Cardirex 3T/4/6T/62,Cathcor, Megacard ElemaA/B,R,R/E,Mingocard 3/4/7,Mingograph 4/62/81 Mingograph 82/740,Mingolog7,Mingophon3/7/310/410/420/710/720,Miracard 1, Sirecard440/440S/460/460S/740Siereg, Sirecard 803,Sirecust 701; SICARD 440
EKG187AOD		integrated; AHA,DB15M-Din3.0		
EKG188IOB		integrated; IEC,DB15M-Banana4.0		
EKG189AOB		integrated; AHA,DB15M-Banana4.0		
EKG190IOS		integrated; IEC, DB15M-snap		
EKG191AOS		integrated; AHA, DB15M-snap		
EKG192IOC		integrated; IEC, DB15M-clip		
EKG193AOC		integrated; AHA, DB15M-clip		

### Medtronic/Physio Control

Part No	OEM No	Description	Picture	Compatible
EKG194IOD		integrated; IEC, 12 pin-Din3.0		Lifepak 12
EKG195AOD		integrated; AHA,12 pin-Din3.0		
EKG196IOB		integrated; IEC,12 pin-Banana4.0		
EKG197AOB		integrated; AHA,12 pin-Banana4.0		
EKG198IOS		integrated; IEC, 12 pin-snap		
EKG199AOS		integrated; AHA, 12 pin-snap		
EKG200IOC		integrated; IEC, 12 pin-clip		
EKG201AOC		integrated; AHA, 12 pin-clip		

### Dixtal EKG cable


Part No	OEM No	Description	Picture	Compatible
---------	--------	-------------	---------	------------

EKG202IOD		integrated; IEC, DB15F-Din3.0		<i>EP-3; DX-203</i>
EKG203AOD		integrated; AHA,DB15F-Din3.0		
EKG204IOB		integrated; IEC,DB15F-Banana4.0		
EKG205AOB		integrated; AHA,DB15F-Banana4.0		
EKG206IOS		integrated; IEC, DB15F-snap		
EKG207AOS		integrated; AHA, DB15F-snap		
EKG208IOC		integrated; IEC, DB15F-clip		
EKG209AOC		integrated; AHA, DB15F-clip		

### Welch Allyn EKG cable

Part No	OEM No	Description	Picture	Compatible
EKG210IOD		integrated; IEC, DB15M-Din3.0		<i>Welch Allyn Cardioperfect MD</i>
EKG211AOD		integrated; AHA,DB15M-Din3.0		
EKG212IOB		integrated; IEC,DB15M-Banana4.0		
EKG213AOB		integrated; AHA,DB15M-Banana4.0		
EKG214IOS		integrated; IEC, DB15M-snap		
EKG215AOS		integrated; AHA, DB15M-snap		
EKG216IOC		integrated; IEC, DB15M-clip		
EKG217AOC		integrated; AHA, DB15M-clip		
EKG218IOD		integrated; IEC, Lemo 14P-Din3.0		<i>Welch Allyn Cardioperfect Pro</i>
EKG219AOD		integrated; AHA,Lemo 14P-Din3.0		
EKG220IOB		integrated; IEC,Lemo 14P-Banana4.0		
EKG221AOB		integrated; AHA,Lemo 14P-Banana4.0		
EKG222IOS		integrated; IEC, Lemo 14P-snap		
EKG223AOS		integrated; AHA, Lemo 14P-snap		
EKG224IOC		integrated; IEC, Lemo 14P-clip		
EKG225IOC		integrated; AHA, Lemo 14P-clip		

## Nihon Kohden

Part No	OEM No	Description	Connector	Compatibility
EKG226IOD		integrated; IEC, DB15M-Din3.0		
EKG227AOD		integrated; AHA, DB15M-Din3.0		
EKG228IOB		integrated; IEC, DB15M-Banana4.0		
EKG229AOB		integrated; AHA, DB15M-Banana4.0		
EKG230IOC		integrated; IEC, 10 leads clip end		
EKG231AOC		integrated; AHA, DB15M-clip		
EKG232IOS		integrated; IEC, DB15, snap		
EKG233AOS		integrated; AHA, DB15M-snap		

## Philips ECG cable, new ten leads

Part No	OEM No	Description	Connector	Compatibility
EKG234IOD		integrated; IEC, DB15M-Din3.0		
EKG235AOD		integrated; AHA, DB15M-Din3.0		
EKG236IOB		integrated; IEC, DB15M-Banana4.0		
EKG237AOB		integrated; AHA, DB15M-Banana4.0		
EKG238IOC		integrated; IEC, 10 leads clip end		
EKG239AOC		integrated; AHA, DB15M-clip		
EKG240IOS		integrated; IEC, DB15, snap		
EKG-10026H		integrated; AHA, DB15M-snap		



EKG-100261		ECG trunk cable ten leads		
------------	--	---------------------------	--	--